

M155

Intramuscular Injection model of Upper arms



SAKAMOTO MODEL
CORPORATION

User's Manual (151218)

— I N D E X —

- Introduction/Table of Contents 1
- Purpose/Content 2
- Precautions concerning handling 3
- Usage/Storage method 4

! Note	Please open your parcel and check the main unit and accessories of the product immediately upon receipt. If the products are delivered damaged or with parts missing, please contact our company or the distributor.
-------------------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

株式会社坂本モデル



Introduction

Thank you for purchasing our product, “Intramuscular Injection simulator of Upper arms”

Read this user’s manual carefully in order to use the product for a long time appropriately.

This product is closely checked before factory shipment. However, when this product arrives, open it soon, and check if anything on the main unit and accessories are missing or broken.

Table of contents

Introduction / Table of Contents.....	1
Purpose / Content.....	2
Precautions concerning handling.....	3
Usage / Storage method.....	4

Purpose

Training can be performed for the following items:

1. Judgment of intramuscular injection of Upper arms

Contents

1. Main unit (including the drainage bag)
2. Battery cells (four AA alkali cells) 1 set



Precautions concerning handling

1. When performing practice using this product, use a fresh needle tip. Otherwise, the epidermis and/or sensor might be degraded.
2. Use water (tap water) only for training. Do not use any fluids other than water.
3. Pay sufficient attention to handling the product because it is fragile.
4. Do not apply a load on the transparent section. Doing so might damage the product.
5. When the model gets dirty, wipe off the dirt using a soft cloth with a neutral detergent. (Do not use benzine, thinners, or polishing powders.)
6. Since the cover might drop if raising the model by holding it on the cervix, hold the model under both arms to raise it.
7. Use four AA alkali battery cells.
Set them while paying attention to polarity. When the volume of the buzzer becomes smaller or the lamp becomes darker than before, replace them with new ones.
8. If anything strange is detected when this model is being used, stop using it and notify the dealer or us as soon as possible.
9. This product is a simulator for users to perform advanced training. Nonetheless, it does not guarantee that a similar effect can be obtained by actual treatment.
10. Note in advance that this product is subject to change without notice to improve it.

Usage

- 1) Set four AA alkali battery cells in the battery box while paying attention to polarity.
- 2) Turn on the power switch.
- 3) The injection site can be confirmed by palpating the acromion.
- 4) Perform training.
- 5) It is judged by the embedded lamp and buzzer whether intramuscular injection has succeeded or not.
 - If training was performed in a proper position including depth, the green lamp will light up.
 - If the injected site was not appropriate, the red lamp will light up, and the buzzer will sound.
 - If that position was correct but the needle was inserted too deep, both the red and green lamps will light up, and the buzzer will sound.

Storage method

1. After the fluid is drained, stand the model and store it.
2. Most of the water injected in an appropriate position will be drained after a while.
Water drops left in the injected sites or the drain tube will not leak out but dry soon because the material allows air to pass through.
3. Do not store the product in a place exposed to direct sunlight or hot and humid air.
4. Be sure to take out the battery cells before storing the model.
Be careful as the battery cells might corrode and cause a power fault if not removed.

SAKAMOTO MODEL
C O R P O R A T I O N

Head office

34,Higashitakagicho,Shimogamo
Sakyo-ku, KYOTO 606-0865, JAPAN

Phone: +81-75-701-1135

Fax: +81-75-722-5638

E-mail:info@sakamoto-model.com

URL:http://sakamoto-model.co.jp

Overseas Sales Department

1-9-29,Senbahigashi,Minoh City,Osaka
562-0035, JAPAN

Phone: +81-72-729-8350

Fax: +81-72-729-8385

E-mail:trade@sakamoto-model.com