

Aged Simulation Suit 2

NEW

M176-6 Aged Simulation Suit II

S size

M176-7 Aged Simulation Suit II

M size

M176-8 Aged Simulation Suit II

L size

M176-9 Aged Simulation Suit II

LL size

● Weight/About 7kg ● Accessories/Exclusive goggles(1 pair), Ear plugs(40 pieces), Gloves(10 pairs), Weights ankle(4 pieces), Weights wrist(4 pieces), Cane(1 piece), Belt cover(100 pieces) Storage bag(1 piece)

Easy to put on and take off, ensuring smooth practice

As we age, our bodies undergo various changes that greatly affect our daily lives.

"Elderly Experience Suit II" is a hands-on educational material that simulates these physical changes associated with aging. It is a suit type that is easy to put on and take off, making it ideal for "group trial training" where multiple people can experience it. By experiencing physical changes, you can understand the feelings of the elderly, and please use this experience for care training and product development.



point!

Step-limiting belt added!

《Renewal point!》

Characteristics of the elderly, such as step limiting, knee flexion and stooping, can be simulated. The movement of an elderly person's knee can be more accurately reproduced.



point!

Renewal of goggles!


《Renewal point!》

The goggles have been changed to a glass type, to eliminate the discomfort which many people felt during use.




Features

1 Fit-on



It takes only about 3 minutes to put on the suit. This body suit type gear worn like overalls allows students to experience restrictions in physical functions with little loss of time.

2 Adjustment



- ① The angle of hip flexion is adjusted.
- ② The motion ranges of parts of the body can be adjusted.

The movable range and bendability of the hips can be freely adjusted. Characteristics of the elderly such as stooped backs can be simulated by changing the balance between the lengths of the left and right belts.

3 Restriction Belt



Neck restriction belt


Elbow restriction belt

Knee restriction belt

Newly developed neck/elbow/knee control belts restrict motion ranges.

The wearer can experience the fear of not being able to promptly cope with danger.

The product designed on the basis of data allows the wearer to experience physical changes due to aging



Vision



Hearing



Grip strength



Physical abilities

Specifications

S size	Recommended stature	145cm~155cm
	Weights	200g arm weights (4 pieces), 400g leg weights (4 pieces)
	Body color	Green
M size	Recommended stature	155cm-165cm
	Weights	200g arm weights (4 pieces), 400g leg weights (4 pieces)
	Body color	Pink
L size	Recommended stature	165cm-175cm
	Weights	250g arm weights (4 pieces), 500g leg weights (4 pieces)
	Body color	Blue
LL size	Recommended stature	175cm-185cm
	Weights	300g arm weights (4 pieces), 600g leg weights (4 pieces)
	Body color	Green

Exchangeable parts

Ear plugs 20 pieces
Gloves 10 pairs
Belt cover 100 pieces
special shoes
Special goggles
Cane