


"New TOMOKO" - Nursing Training Manikin

● Accessories/Wig for shampoo, Pajamas, Lubricant oil (Silicone oil) , Baby powder, Catheters, "Veinmate" 2

Good flexibility and soft texture facilitate realistic nursing practice.

Features

- Posture change
- Intramuscular injection
- Intravenous injection /blood sampling
- Oral hygiene
- Gastrostomy care
- Total parenteral nutrition
- Stoma care
- Tracheotomy care
- Changing pajamas
- Transportation
- Bed bath
- Suctioning
- Catheterization
- Nutrition (oral/nasal)
- Enema
- Manual disimpaction
- Intestine lavage/Gastric lavage
- Assisting excretion
- Restriction
- Fomentation
- Bandaging
- Shampoo



Specifications

| | |
|----------|------------|
| Height | 160cm |
| Weight | About 17kg |
| Material | Silicone |

Accessories

- "Veinmate" 2 (one unit only including a Irrigator and a Belt)
 - Wig for shampoo
 - Pajamas
 - Lubricant oil(Silicone oil)
 - Baby powder
 - Catheters
 - Balloon catheter 14Ft
 - Parenteral catheter No. 4
 - Nelaton catheter 8EG
 - Gastric catheter 10mm
 - Nelaton catheter 15EG
- *Training equipment other than the above five items is not included in this product.



! "Veinmate" 2 (one unit) is included in New TOMOKO.

A wide range of practices based on the nursing education program is possible.

New Tomoko is designed based on the nursing education program, using in various practical subjects is possible. Trainees can acquire a wide range of skills essential in clinical practice, such as the basic nursing assistant skills to patient's daily living needs as well as care skills associated with catheter insertion.

Practice

Tracheotomy care



Suctioning



Gastrostomy care



Total parenteral nutrition



Nutrition (oral/nasal)



Intestine lavage/gastric lavage



Urinary catheterization



Enema administration



Stoma care



Posture change



Oral hygiene



Intravenous injection



*Intravenous injection practice can be practice by using included "Veinmate" 2.

M102-1 Interchangeable to male genitalia model

Male Genitalia Model for "New TOMOKO" (M102)

- Weight/0.5kg/unit
 - Accessories/Irrigator, Intestines
- <Irrigator stand sold as a separate item>

