


Sakamoto Acupuncture Simulator for Sciatic Nerve

● Weight/About 3.2kg

Practical training using a lateral position model simulating an actual patient.

This model permits the practical training at a lateral position simulating actual treatment, aiming at mastering the accurate acupuncture skills for each individual postures. Since the skin is made of special material allowing more practical training close to the actual environment. The iliac crest (iliac spine), great trochanter, and ischial tuberosity are incorporated in the model. The trainees can practice palpating, and the green light turns on and a buzzer sounds when the needle is inserted to the correct point facilitating easy self-practice.



Practice

1. Technique 1




2. Technique 2



3. Needle insertion



Posterior superior iliac spine - great trochanter



Ischial tuberosity - great trochanter

White-colored Educational Manikin

● Size/H160cm

Educational model of the human body on which writing with a water-soluble ink is possible.

