

# M212-1/M212-2

## Transparent Female Catheterization & Enema Model Transparent Male Catheterization & Enema Model

●Weight/approx. 8.5kg ●Case size/W58×D46×H29cm

●Accessories/Bottle(1L), Bottle stand, Tray, Balloon catheter 16Fr, Nelaton catheter 8EG, Lubricant oil, Baby powder, Storage case

Bladder and rectum are transparent, allowing confirmation of catheter insertion status. The prostate is also realistically reproduced.



Urination training enables actual outflow of liquid (water).



Urination training enables actual outflow of liquid (water).

### Feature

●External genitalia closely resemble real anatomy. When the labia minora are opened, mucosal folds and urethral opening can be observed.



●Balloon catheter and enema insertion techniques can be practiced, providing realistic resistance and pressure sensations. The prostate and rectum are newly reproduced.



Note: External genitalia cannot be removed from the main model.

●Pelvis is included, allowing practice of catheterization angles based on bladder position.

