

M148-3

"Veinmate" 2- I.V. Injection Pad

- Quantity/5 pieces/set
- Weight/About 0.7kg (5 pieces)
- Case size/W29×D20×H10cm
- Accessories/Irrigator, Belt, Storage case


Practices of different levels

- Makes possible practice closely simulating clinical situations in various patients.
- The epidermis is very close to that of a living body and leaves no needle mark even after countless practices.




Feature

The feeling of the skin and vessels on palpation and needle insertion has been made closer to that in the human body.



- The vessels are arranged in meandering courses.



- The vessels pass different depths from close to the surface to deep areas.

Practice

Preparation



Palpation



Intravenous injection



Blood sampling



M117-8

Stock Solution for Simulated blood 100 cc

M148-2

Irrigator

- Y connector included