

# M161

## Resusci Anne Skill Reporter Model

- Case size/W90×D32×H51cm
- Accessories/Skill Reporter, Training mat, Pulse Bulb with Tubing, Airways(2 pieces), Adult Manikin Face Skin(3 pieces), Manikin Wipes(50 pieces), Storage case

With rich support functions, the model can also be used for CPR practice by medical workers.

- The 3-color representation of the bar graph reflects the intensity level of artificial respiration and cardiac massage.
- Easy-to-perceive displays using illustrations of the human body.
- The cord connecting the model to the Skill Reporter has a sufficient length.
- Automatic event recording by the Skill Reporter  
The sensors built in the model automatically records the following items.  
①Checking the level of consciousness,  
②Securing of the air way, ③Feeling for pulse,  
④Too rapid blowing in, ⑤Incorrect hand position



(Skill Reporter)

# M162

## Resusci Anne Skill Guide Model

- Case size/W90×D32×H51cm
- Accessories/Skill Guide, Training mat, Pulse Bulb with Tubing, Airways(2 pieces), Adult Manikin Face Skin(3 pieces), Manikin Wipes(50 pieces), Storage case

With easy-to-understand CPR evaluation functions

- Determination of the correct site of cardiac massage.
- The Skill Guide has indicators of the tidal volume, intensity of chest compression, position of chest compression, and stomach distension and checks whether CPR training has been conducted appropriately.



(Skill Guide)

# M122-6

## Resusci Junior

- Case size/W73×D26×H44cm
- Accessories/Skill Guide, Training mat, Pulse Bulb with Tubing, Airways(3 pieces), Junior Manikin Face (3 pieces), Manikin Wipes(18 pieces), Storage case

Model of a 5-year-old child to effectively learn resuscitation techniques on children.



(Skill Guide)