


# Transparent Suction Model

- Weight/About 5.5kg ● Case size/W68×D34×H32cm
- Accessories/Suction tube with LED light, Lubricant oil, Storage case

It is “Transparent” to understand anatomical structure. Oral suction, nasal suction, tracheostomy cannula suction possible.



## Practice

### 1 Lifting up the head part



### 2 Moving the head part to left and right



## Feature

## 1 Can observe the inside.

You can observe tongue, epiglottis, trachea and esophagus as it is transparent from head to trachea. You can learn the nasal/oral structure and relation between trachea and esophagus.



▲Exclusive suction tube

## 2 Can check the tube location by using catheter with LED light.

Catheter with LED light comes with model. You can observe the catheter moving when suction.



Suction position of pyriform sinus



Suction position of bronchi

## 3 Can check suction tube at various location.



Practice using pseudo sputa is not possible.