

M175

Sakamoto Suction Simulator


- Weight/About 15kg ● Case size/W68×D34×H32cm
- Accessories/Pseudo sputa, Baby powder, Lubricant oil, Container for mixing pseudo sputa, Muddler for mixing pseudo sputa, Syringe with tube, Silicone cap, Flushing funnel, Flushing tube, Tray

Basis of respiratory management
Endotracheal intubation can be learned.



Feature

The head can be moved to the left and right



This model is a simulator with a lifelike texture, and various suction methods can be learned.

The suction of secretions not only in the oral/nasal cavities but also in tubes inserted via the oral cavities is possible.

Since the trachea has already been incised, suction with tracheal cannula insertion can also be performed.

Various types of suction training are possible

Suction of secretions in the oral cavity



Suction of secretions in inserted tubes



Suction of secretions in the nasal cavity



Suction of secretions in the tracheal cannula



Washing

Washing of the bronchi



Washing of the oral cavity



Exchangeable parts

Pseudo sputa