

# M179

## Vital Signs Simulator Sakamoto Baby Touch

- Weight/7.2kg (Including accessories) ● Case size/W55×D32×H37cm
- Accessories/Baby cloth, Stethoscope, Thermometer, Baby powder

Notice the symptoms promptly!  
Represent the neonatal vital signs precisely.




### Feature

- 1 Available in four languages!  
(Jap, Eng, Chin, Kor)



Language select screen



English



Chinese




Korean

- 2 Can set vital signs (Heart rate, Pulse rate, Respiration rate, Body temperature) at 5~6 levels.

- 3 Simply set with touch panel



- 4 You can practice measuring the pulse of anterior fontanel, and checking of cryptorchidism (undescended testicle)




Renewal



## Heart Rate

Confirming the heart sound near the heart

Can be set at 6 levels:  
"60·80·120·140·160·200"



## Pulse Rate

Measuring the pulse at 6 points

Can be set at 6 levels:  
"60·80·120·140·160·200"



Superficial temporal arteries



Carotid arteries



Radial arteries



## Respiration Rate

By up and down movement of abdominal area

Can be set at 6 levels: "20·40·50·60·80·100"



## Body Temperature

Rectum

Can be at 5 levels: (About 32°C - 42°C)

