

M181

Fit-on Intravenous Injection Trainer I.V.Pad

- Quantity/5 pieces/set ● Weight/About 1.3kg (5 pieces) ● Case size/W45×D27×H10cm
- Accessories/Exchangeable blood vessel (normal) 10 pieces, Stick for exchange vessels 1 piece, Storage case

Easier and speedier training of injection and blood sampling!



Feature

Just two seconds to fit-on!

1 Open!




2 Wrap tightly!



Easy vessel change. Easy maintenance.

1 Connect



2 Pull out!



Practice

Fits better to arm



Palpation




Intravenous Injection

Insertion of indwelling needle
Fixing with adhesive tape

Exchangeable parts

Two types of training is possible using blood vessels sold separately.



Exchangeable blood vessels (normal) 100pcs
Exchangeable blood vessels (narrow) 100pcs

< Stock solution for simulated blood 100cc (sold separately) >

M182

Arm for Simulator

● Size/W60×D25×H23cm

*M182 can be used for M118-1, M148-3, M178, M181, M183, AM-1, AM-2

