


"New TOMOKO" - Nursing Training Manikin

- Accessories/Washable wig, Pajamas, Lubricant oil, Baby powder, Catheter set, I.V.Pad, Simulated gastrostomy catheter, Simulated tracheal cannula, Exclusive tape (25mm width x 5m, 50mm width x 5m) 1 each

Good flexibility and soft texture
facilitate realistic nursing practice

Features

- Posture change
- Intramuscular injection
- Intravenous injection /Blood sampling
- Nutrition (oral/nasal)
- Gastrostomy care
- Central Venous nutrition
- Stoma care
- Tracheostomy care
- Changing pajamas
- Transportation
- Bed bath
- Suctioning
- Urinary catheterization
- Nutrition (oral/nasal)
- Enema
- Manual Disimpaction
- Intestine lavage/Gastric lavage
- Assisting excretion
- Restriction
- Fomentation
- Bandaging
- Shampoo



point!

RENEWAL POINT

The list of things that can be practiced have increased.



Tracheal cannula



Exclusive tape



Dedicated tape



Immobilization of joints

Specifications

Height	165cm
Weight	About 17kg
Material	Silicone

Accessories

I.V.Pad (1 piece only)

- Accessories/Blood bottle

Exclusive tape (25mm width x 5m, 50mm width x 5m)

1 piece each -Total 2 sets -

Washable wig

Pajamas

Lubricant oil

Baby powder

Simulated gastric fistula catheter

Simulated tracheal cannula

Catheter set

- Balloon catheter 14Fr
- Gastric catheter 10mm
- Feeding catheter No. 4
- Nelaton catheter 15EG
- Nelaton catheter 8EG



! I.V.Pad (1 piece) is included in New TOMOKO



Simulated gastric fistula catheter



Exclusive tape (2type)

A wide range of practices based on the nursing education program is possible

New Tomoko is designed based on the nursing education program, being used in various practical subjects is possible. Trainees can acquire a wide range of skills essential in clinical practice, such as the basic nursing assistant skills to patient's daily living needs as well as care skills associated with catheter insertion.

Practice

Tracheotomy care



Suctioning



Gastrostomy care



Central venous nutrition



Nutrition (oral/nasal)



Intestine lavage/Gastric lavage



Urinary catheterization



Enema administration



Stoma care



Posture change



Oral hygiene



Intravenous injection



* Intravenous injection practice can be practice by using included I.V.Pad

M102-1

Interchangeable with male genitalia model

Male Genitalia Model for "New TOMOKO"



- Weight/1 piece 0.5kg
 - Accessories/Irrigator, Intestines
- <Irrigator stand is sold separately>

M101-5

Storage bag for nursing doll for one doll

- Material/Nylon

