

Wearable-type Wound Trainer

- Accessories/ Transparent attachment belt, Exclusive stand
- Unit weight: 0.45kg
- Unit size: W37×D26×H2cm
- Case size: W39×D28×H5cm
- Material: Soft urethane

6 types of reusable wounds are available to choose from!



Features

Sakamoto Model's New Approach!!

They can be attached to the human body or onto a mannequin. Plus the viscosity returns simply by being washed in water, so that repeated use for training is also possible.



Practice

1 Reproduce wounds on humans and mannequins in many different places



Human



Mannequin



2 Antiseptic care while simulating actual wounds



3 Suitable for using with other nursing techniques (when replacing nightwear)



point!

For places which are difficult to attach



The special accessory belt can be used to allow practising in places where it is difficult to attach to actual skin surfaces.