

M208

NEW

Sakamoto Intravesical Ultrasound Observation Trainer

- Weight/about 1.2kg ●Size/H11.5×W12×D13cm
- Accessories/Lubricant oil, Injection syringe, Baby powder, Cup (300cc), Exclusive balloon catheter 16Fr, 2 Trays

Observe in real time, changes in urine residual volume as you urinate



Ultrasound image during balloon placement

Ultra sound

M210

NEW

Radial Artery Line Securing Trainer (Ultrasound compatible)

- Case size/W28×D32×H23cm
- Accessories/10 Exchangeable blood vessels, Stick for exchange vessels, Cup

Securing an arterial line
under echocardiography -
Compatible with training



Puncture under ultrasound



Exchangeable parts

- M210-1 Exchangeable blood vessel for radial artery line securing trainer
- M210-2 Exchangeable pad for radial artery line securing trainer