

Female Catheterization & Enema Simulator Male Catheterization & Enema Simulator

- Weight/About 5.5kg ● Case size/W62×D44×H30cm
- Accessories/Simulated soft stools (yellowish-brown), Simulated soft stool mixing container, Simulated soft stool stirring rod, Tray, Replacement non-return valve (×2), Bottle, Bottle stand, Balloon catheter 16Fr, Nelaton catheter 8EG, Lubricant oil, Baby powder, Storage case, Foreskin (men only)

Catheterization and enema practice can be carried out using a single model



1 Removable



Genital parts can be easily removed at a single touch. The main body can be used for training by switching the male and female genitalia.

2 Improved durability



Updated materials for genitals which could previously break easily. Texture is closer to that of a real body, while durability has been improved.

3 Possible to place balloon catheter



Catheters from 12Fr to 18Fr can be used.
The state of balloon catheter insertion can be checked.

4 Simple design so that the bottle can be easily removed and emptied after use



M200-3

Catheterization & Enema Simulator (Male Female genitalia set)

- Configuration/1 body, 1 male genitalia, 1 female genitalia
- Accessories/Simulated soft stools (yellowish-brown), Simulated soft stool mixing container, Simulated soft stool stirring rod, Tray, Replacement non-return valve (×2), Bottle, Bottle stand, Balloon catheter 16Fr, Nelaton catheter 8EG, Lubricant oil, Baby powder, Storage case, Foreskin

