

# M169

## Fit-on Manual Disimpaction Simulator

- Weight/About 4kg ● Case size/W55×D26×H27cm
- Accessories/Belt, Simulated excrement, Lubricant oil, Baby powder, Storage case
- Director/Hiromi Tashiro

Exclation



### Practice

Set Simulated excrement can be easily set and used repeatedly



1.Remove the cap on the top of the main unit. It is easy to pull out by pinching the dent.



2.Take a moderate amount of simulated stool and put it inside the body (rectal area) Insert it by pushing it inside.



3. Close the cap so that it is aligned with the surface of the main unit. Please close the lock at the top.

Easy to install Easy to install by simply fixing it with a belt and fitting it on the body



Since it is wearable, you can practice in the actual used positions. Suitable for pair or group learning.

The anal area can be removed for cleaning



### Exchangeable parts

Simulated excrement  
Lubricant oil