

Pregnancy Experience Suit

NEW

- Weight/About 8.7kg, Abdominal weight: About. 6kg (750g × 8 pcs),
Breast weight: about. 1.5kg (750g × 2 pcs)
- Material/Suit: Nylon 100%; Weight bag: Nylon 100%;
Weight contents: Iron, urethane
- Accessories/Weights (750 g × 10 pcs), Storage bag × 1
<Apron sold separately>

Experience 24 hours a day what it is like to be pregnant with the Pregnancy Experience Suite

The changes that occur during pregnancy are not limited to just a swollen abdomen. As well as the abdomen growing bigger with the new life forming inside it, there are of course other changes occurring both physically and mentally, and they have an effect on all aspects of daily life. In order for everyone to experience the changes involved with pregnancy. In order to make it possible to practice in a clean state. The lifelike design of Sakamoto Model's Pregnancy Experience Suit is the ideal product for you to use.



Features

1 Free design size which can be fitted easily to either men's or women's bodies



Cord lengths can be adjusted so that the suit can be fitted easily to either men's or women's body styles. The simple design lets the suit fit itself just by tightening the buckle.

3 Buckle type which does not get caught in hair or clothing



Uses Sakamoto's original fitting method. This eliminates any discomfort from the product's belt fittings getting caught in hair or clothing.

2 Adjustable weights



The weights in the breasts and abdomen can be inserted and removed so that you can adjust the total weight during training.

4 Washable so it can stay clean at all times



The experience suits may be worn by many different people. Because the suit itself can be easily washed, it can be kept clean so that you can always use it with confidence. (Remove the weights from the suit before washing.)