

"Musclemate" 2 - Intramuscular Injection Simulator

● Weight/About 0.5kg ● Case size/W32×D12×H11cm ● Accessories/Drain bag, Belt, Storage case

Mastery of the hard-to-learn feeling of "depth"



Features

Correct or incorrect is indicated with built-in lamps and a buzzer.



[Green lamp]

OK

The Green lamp is lit when the injection site and the depth are correct.



[Red and green lamps]

Error regarding the depth

Both the red and green lamps are flashed, and the buzzer goes on, when the injection site is correct but too deep.



[Red lamp]

Error regarding the site

Both the red and green lamps are flashed, and the buzzer goes on, when the injection site is correct but too deep.

* The light of the LED lamps can be seen through the skin.

Preparation



Easy installation with Velcro tape.

Palpation



A built-in acromion ensures proper confirmation of the injection site.

Pinching



When you pinch it with your fingers, the muscles swell realistically and you can determine the injection site.

Brachial I.M. injection



It is possible to measure the correct injection site including the depth. The epidermis is extremely close to the real body and needle marks do not remain even after repeated practice.

Subcutaneous injection



Subcutaneous injection is also possible.

M182

Arm for Simulator

● Size/W60×D25×H23cm

※M182 can be used for M118-1, M118-2, M148-3, M178, M181, M181-3, AM-1, AM-2

