

"Musclemate" Air (5 pieces per set)

NEW

- Weight/About 150g ●Size/W10×D30×H7cm
- Components/"Musclemate" Air main unit (5 pieces), 5 attachment belts (5 pieces)

Light, durable, and high quality! We offer a 5 piece set at this price!!
Ideal for a series of intramuscular injection operations such as vaccinations.

- With a structure that does not use electricity, it is ideal for training that assumes the administration of physiological saline and various drugs.
- The epidermis is made of a special material that does not leave needle marks, so it can be used for repeated insertions and is replaceable.
- The internal structure reproduces the acromion, making it ideal for procedures to confirm the correct insertion site.
- The mounting surface is made of a hard material that does not allow injection needles to pass through, and is designed to fit people of various body types. It can be secured with a special belt.



Features

Ideal for large group vaccination training!!



Injection pad can be easily removed.



The acromion is also built-in, allowing you to check your technique.



Drainage is easy, and chemical solutions can be used as well.