

"Veinmate" 2- I.V. Injection Pad

● Quantity/5 pieces/set ● Weight/About 0.7kg (5 pieces) ● Case size/W28×D22×H11cm ● Accessories/Irrigator, Belt, Storage case

Realize 4 types of training by arranging “thick/thin” blood vessels and “shallow/deep” points


- Makes practice possible closely simulating clinical situations on various patients
- The epidermis is very close to that of a living body and leaves no needle mark even after countless practice.

Features

The feeling of the skin and vessels on palpation and needle insertion has been made closer to that of the human body



- The vessels are arranged in meandering courses.



- The vessels pass different depths from close to the surface to deep areas.

Practice

Preparation



Easy installation with Velcro.

Palpation



Reproduces both the parts where the blood vessel stands out and the part of the blood vessel deep in skin. You can check the blood vessel position by touching it with your fingertip.

Intravenous injection



You can practice on the skin and blood vessels with realistic feeling.

Blood sampling



There are two types of blood vessels, thick and thin which can be used according to the level of the practice.

M117-8

Stock solution for Simulated blood 100cc

M148-2

Irrigator

- Y connector included

M182

Arm for Simulator

- Size/W60×D25×H23cm

※M182 can be used for M118-1, M118-2, M148-3, M178, M181, M181-3, AM-1, AM-2

