



"Buttockmate" type 2

● Weight/About 1.2kg ● Case size/W45×D15×H26cm ● Accessories/Belt, Drain bag, Storage case

Reproduces Clarke's point and Hochstetter's point

A wider range of practical training is now possible while maintaining the greatness of the "Buttockmate". Practicing can be done in a pair of two people, and it is possible to take the actual position (lateral recumbent position).

Because it is a wearable type, the patient's anxiety can be felt, and it can also encourage growth in students in terms of "caring".



Selection of an incorrect site is indicated by the blinking of the red lamp and the sound of buzzing. Successful puncturing at the correct site is indicated by the green lamp.
(Only at Clarke's point and Hochstetter's point)



The epidermis is made of a special material that leaves no needle marks.

Being compact, the device needs only a small space for storage.

Practice

Hochstetter's point



Skeletal reference points such as the anterior superior, anterior posterior iliac spines and greater trochanter can be felt for efficient training in finding the injection site.

Clarke's point



Fluid can be injected.

* Do not use a drug solution. Use water (tap water only) for practice.