

“Buttockmate” Right

Intramuscular Injection Simulator

NEW

●Weight/About 1.8kg ●Case size/W46×D29×H16cm ●Accessories/Belt, Exclusive stand, Storage case

Practice on Clark's point, Hochstetter's site Upper outer quadrant of buttocks

Supports three typical injection methods. The soft skin unique to the Sakamoto model and the internal sensors are all replaceable. High-level trainings can be realized at a low running cost.



Also possible to practice by using exclusive stand.

Exchangeable parts

Skin cover
(Absorbing) Sponge
Sensor
Wearable belt

Features

The parts are consumable and can be replaced




The absorbing sponge, sensor and skin cover can all be removed and replaced. This product is easy to maintain.

Soft and realistic skin



The feel of both the skin and skeleton are faithfully reproduced. Of course it is also possible to use tap water for injection.

Judgement indicator when successful



When the injection is carried out correctly, the judgement indicator will light up green and a audible chime can be heard.