

"Buttockmate" - Intramuscular Injection Simulator

● Weight/About 2kg ● Case size/W46×D29×H16cm ● Accessories/Belt, Drain bag, Transparent package, Storage case

Practice gluteus injection skill by palpation



The epidermis is made of a special material that leaves no needle marks.



A drawing of the injection site is written on the transparent package.

- All students can participate in the class at the same time. Students can practice many times during one class contributing to further improvement of their skills.
- Students can train in pairs playing both the role of nurse and that of the patient alternating. This wearable type enables comparison with a real body.
- As ilium crista, greater trochanter and iliac spine are built in, confirmation and practice for the gluteus injection skill by palpation is available.
- The flashing of a red light and the sound of the buzzer warn in case of a missed injection, incorrect injection point or a too deep injection. The green light turns on in case of a correct injection point.
- Because it is a wearable type, you can attach it to the human body or model doll for more realistic and practical practice.

Practice



By palpation, injection training of the correct site is possible.



Water (only tap water) can be injected.



You can also practice by yourself using a package.