

Fit-on Intravenous Injection Trainer I.V.Pad

- Quantity/5 pieces/set ● Weight/About 1.3kg (5 pieces) ● Case size/W45×D27×H10cm
- Accessories/Exchangeable blood vessel (normal) 10 pieces, Stick for exchange vessels 1 piece, Storage case

Easier and speedier training of injection and blood sampling!



Features

Just two seconds to fit-on!

1 Open!



2 Wrap tightly!



Easy vessel change. Easy maintenance.

1 Connect



2 Pull out!



Practice

Fits better to arm

By expanding the injection range that comes into close contact with the skin, both the injector and the patient can feel even more fear.



Palpation



Intravenous Injection

Insertion of indwelling needle
Fixing with adhesive tape

It is also possible to create meandering blood vessels within the epidermis. You can create various cases.

Exchangeable parts

Two types of training is possible using blood vessels sold separately

You can increase the difficulty by using small blood vessels.

<Normal>



4mm

<Narrow>



3mm

M181-1 Exchangeable blood vessels (normal) 100pcs

M181-2 Exchangeable blood vessels (narrow) 100pcs

<Stock solution for simulated blood 100cc (sold separately)>

M182

Arm for Simulator

● Size/W60×D25×H23cm

※M182 can be used for M118-1, M118-2, M148-3, M178, M181, M181-3, AM-1, AM-2

