

# Sakamoto Electric Venipuncture Trainer

NEW

Injection

- Weight/ 4.5kg ● Case size/W61×D28×H20cm
- Accessories/Arm model ×2, Injection pad (with blood vessels) ×2, Tank with integrated pump, AC adapter, Exchangeable blood vessels ×10 (standard size), Stick for exchange vessels, Baby powder, Spare stopper cap ×2, Cup

A new double lineup for training in intravenous injection is born  
Low cost and easy to use with solutions



## Features

### Low-cost, high-quality injection pads

This product is softer than previous pads and seeks to provide a more realistic feel. It is also possible to create backflows and to use vacuum blood collection tubes. (for both M203 and M204)



### Blood flow rate can be freely adjusted

The flow rate can be adjusted according to the number of arms used and the patient's blood flow.



You can freely adjust the flow rate!

### We have developed a tank with integrated pump which is compact and easy to carry

Compact tank-integrated pump makes preparation easy. The water fills the blood vessels and preparation is complete simply by turning on the power for the pump kit, and cleaning up is completed by tilting the arm so that the water flows back into the tank. Either an AC power supply or batteries can be used! You can use it without having to worry even in places where there is no outlet.

## Exchangeable parts

- M203-1 Exchangeable pad for venipuncture trainer
- M181-1 Exchangeable blood vessel (standard)(100 per bag)
- M181-2 Exchangeable blood vessel (narrow) (100 per bag)
- ※These parts can be used for both the M203 and M204.