

M184

Sakamoto Obstetric Assistant Trainer

- Weight/About 13kg ● Case size/W74×D37×H55cm
- Accessories/External genitalia (skin color/transparent), Fetus (Delivery), Lubricant oil, Baby powder, Storage Case

Transparent external genitalia make it possible to discern the status of fetal head



Obstetric/Pregnant

Features

Two features (Dual)

1 Two types of external genitalia (skin color/transparent)

Can see the status of fetal head through the transparent genitalia.



2 Two types of posture

Supine position, all fours possible.

1 Obstetric assist

It is possible to realistically train a series of birth assistance techniques such as protecting the perineum, protecting the baby's head, delivering the scapula anteriorly (delivering the scapula posteriorly), and delivering the trunk. Also, those who are "manipulating the fetus" are required to understand rotation, which will help them learn how to assist.

● Protection of the perineum



● Crowning



● Delivery of anterior shoulder



2 Can check placenta/amnion

The schematic placenta has amniotic membrane attached to it, and a simulation can be performed to check whether it remains in the uterus. Also, the umbilical cord that connects to the placenta is made of soft material and can actually be clipped.



4 Two types of posture

Supine position, all fours possible



▲ All fours (Knee-chest position)

3 Umbilical ligation and umbilical cord cutting

Forceps can now be used as well.



▲ Hemostatis using forceps



▲ Umbilical cord cutting

Exchangeable parts

Exchangeable genitalia (skin color)

Exchangeable genitalia (transparent)

Exchangeable umbilical cord (10 pieces/set)