

M186-1

Sakamoto Whole Body Pregnant Model II

Renewal

● Accessories/ External Genitalia (skin color/transparent), Fetus(Delivery), Lubricant oil, Baby powder, Placenta, Umbilical cord (×1)

Useful for serial training of birth assist!



Specifications

Height	165cm
Weight	about 19kg

Exchangeable parts

Exchangeable genitalia (skin color)
Exchangeable genitalia (transparent)
Exchangeable umbilical cord (10 pieces/set)

Features

- 1 A renewed design, with the whole body arms and legs as well as the abdominal cover being softer



- 2 Can see the status of fetal head through the transparent genitalia



- 3 Can check placenta/amnion

The textured placenta includes an amnion, so that you can simulate its presence or absence inside the womb. In addition, the umbilical cord which is attached to the placenta is made of a flexible material so that it can actually be clipped.



- 4 Umbilical ligation and umbilical cord cutting

Forceps can now be used as well.



Free-style delivery possible




▲Supine position




▲Semi-supine position



▲All fours (Knee-chest position)



▲Squat position



▲Standing position