

## Physical Assessment Doll TOMOKO Plus

- Accessories/Tablet PC, Control box, Exclusive bedside monitor, Stethoscope, Exclusive sphygmomanometer, Exclusive electrodes with leads, Penlight, Exclusive pulse oximeter, Pajamas, Baby powder

### Smart operation    Smooth practice



Group auscultation and practice are now possible!  
It is now possible to output the auscultation sounds of breathing sounds and heart sounds from the tablet.



#### Specifications

Height	160cm
Weight	26kg (Body 15kg, Accessories 11kg)

#### Options

- M103-1 Dedicated stretcher
- M103-2 Exclusive folding bed
- M103-3 dedicated side wagon
- M103-4 dedicated case (for 1 doll)

#### Features


### 1 Starts up quickly just by turning ON the power    Display is clear with simple settings



#### ● Basic menu screen




#### ● Detailed display screen



#### ● Simulation settings can be saved



#### ● Multiple languages available



### 2 Master the basics, covering 8 items in the physical assessment.

Features various symptoms that can be experienced in a medical setting, such as respiration and pupillary reflex. Facilitates training to ensure proper technique and accurate understanding of the patient's conditions during the assessment.

### Check respiratory rate

<Setting> 7 levels ranging from 10-40 breaths/min



### ECG

<Setting> 9 case examples



### Learn pulse palpitation

<Setting> 12 levels ranging from 30-140 beats/min




### Auscultation of cardiac sounds

<Setting> 20 case examples



### Check pupillary reflex

<Setting> 4 types



### point! RENEWAL POINT!

#### Check breath sounds

<Setting> 12 case examples

Auscultation sounds such as heart sounds and breath sounds can be output from the tablet. It is also possible to practice in groups.



### Check pulse oximeter

<Setting> SpO2 80, 90 and 99%



### Blood pressure measurement

<Setting> Range covers 80-200 for the systolic phase (upper) and 50-120 for the diastolic phase (lower)

