

# M178

## Blood Pressure Simulator

- Weight/4.4kg ● Case size/W45×D40×H21cm
- Accessories/ Stethoscope, Storage case

For better measurement technique of pulse palpation and auscultation



※The attached blood pressure monitor is a mercury-free blood pressure monitor. It is also compatible with your existing mercury blood pressure monitor.

### Features

- 1 Attach to the blood pressure and pulse measurement site!
- 2 The pulse palpation on radial artery is possible
- 3 Korotkoff sounds




### Practice

#### Easy setting and Operation

Easy operation with simple and easy-to-understand screen display, buttons, and dials.

Setting Range
Blood pressure/50~250mmHg
Pulse/20~150 beats



It is possible to make indications on monitor invisible. It can be used for examinations.

#### For both pair-training and self-training



It is possible to fit it on "New TOMOKO" (M102) or "Arm for Simulator" (M182).