

Vital Signs Simulator Sakamoto Baby Touch

- Weight/7.2kg (Including accessories) ● Case size/W55×D32×H37cm
- Accessories/Baby clothes, Stethoscope, Thermometer, Baby powder, Storage case

Notice the symptoms promptly!
Represent the neonatal vital signs precisely

Physical



Features

1 Simply set with touch panel



2 You can practice measuring the pulse of anterior fontanel, and checking of cryptorchidism (undescended testicle)



3 Can set vital signs (Heart rate, Pulse rate, Respiration rate, Body temperature) at 5~6 levels

| バイタルサインシミュレーター Baby Touch | | | | Language |
|------------------------------|-----|-----|----|----------|
| 心拍数 | 脈拍数 | 呼吸数 | 体温 | |
| 120 (回/分) | 120 | 40 | 正常 | |
| 60 80 120 140 160 200 | | | | UP |

| バイタルサインシミュレーター Baby Touch | | | | Language |
|------------------------------|-----|-----|----|----------|
| 心拍数 | 脈拍数 | 呼吸数 | 体温 | |
| 120 (回/分) | 120 | 40 | 正常 | |
| 60 80 120 140 160 200 | | | | |

| バイタルサインシミュレーター Baby Touch | | | | Language |
|------------------------------|-----|-----|----|----------|
| 心拍数 | 脈拍数 | 呼吸数 | 体温 | |
| 120 (回/分) | 120 | 40 | 正常 | |
| 20 40 50 60 80 100 | | | | |

| バイタルサインシミュレーター Baby Touch | | | | Language |
|------------------------------|-----|-----|----|----------|
| 心拍数 | 脈拍数 | 呼吸数 | 体温 | |
| 120 (回/分) | 120 | 40 | 正常 | |
| さらに低い | 低い | 正常 | 高い | さらに高い |

4 Available in four languages!

The seamless structure faithfully reproduces the body shape and skin texture of a newborn baby You can practice just like the real thing



Heart Rate

Confirming the heart sound near the heart

Can be set at 6 levels:
"60•80•120•140•160•200"



Pulse Rate

Measuring the pulse at 6 points

Can be set at 6 levels:
"60•80•120•140•160•200"



Superficial temporal arteries



Carotid arteries



Radial arteries



Respiration Rate

By up and down movement of abdominal area

Can be set at 6 levels: "20•40 50•60•80•100"



Body Temperature

Rectum

Can be at 5 levels: (About 32°C - 42°C)

