

M175

Sakamoto Suction Simulator

- Weight/About 15kg ● Case size/W68×D34×H32cm
- Accessories/Pseudo sputa, Baby powder, Lubricant oil, Container for mixing pseudo sputa, Muddler for mixing pseudo sputa, Syringe with tube, Baby powder, Cleaning funnel, Cleaning tube, Silicone cap, Tray, Simulated tracheostomy cannula, Suction catheter 14Fr, Storage case

Suction/Feeding


Basis of respiratory management Endotracheal intubation can be learned



This model is a simulator with a lifelike texture, and various suction methods can be learned.
In addition to intraoral and intranasal suction, oral suction through the intubation tube is also possible.
Since the trachea has already been incised, suction with a tracheal cannula insertion can also be performed.

Features

The head can be moved to the left and right



point!

Pseudo sputa renewal!

The viscosity can be freely adjusted by adjusting the amount of Pseudo sputa.



Various types of suction training are possible

Suction of secretions in the oral cavity



Suction of secretions in inserted tubes



Suction of secretions in the nasal cavity



Suction of secretions in the tracheal cannula



Washing

Washing of the bronchi



Washing of the oral cavity



Exchangeable parts

Pseudo sputa